



PELI
BIO THERMAL™

GOLDEN HOUR PRODUCT RANGE



GOLDEN HOUR™ TECHNOLOGY

Winner of the 2003 "US Army's Greatest Invention Award", Golden Hour™ technology was developed to protect the lives of military personnel on the frontlines by safely transporting blood in a critical amount of time.

The very same patented and award-winning technology provides superior thermal protection for high-value, temperature-sensitive payloads in our renowned Golden Hour™ Mobile and Golden Hour™ Medic containers.

WHY CHOOSE GOLDEN HOUR™ SYSTEMS?

- Ease of use – modular design for quick assembly and efficient inventory storage
- Nimble and tactical for delivery of whole blood to the point of injury within a crucial window of time
- Designed and tested to last in the harshest conditions
- Low profile and life-saving
- Quality, durable construction



GOLDEN HOUR™ MOBILE

DURATION
48 - 72 Hours

SERIES 4 TEMPERATURE
2°C to 8°C

VOLUME
2L

SERIES 22 TEMPERATURE
15°C to 25°C

VOLUME
2L

GOLDEN HOUR™ MEDIC

DURATION
24 - 48 Hour

TEMPERATURE
1°C to 10°C

VOLUME
2L

GOLDEN HOUR™ ONE (F&V)

DURATION
48 Hours


TEMPERATURE
2°C to 8°C

VOLUME
0.5L

GOLDEN HOUR™ PRODUCT SPECIFICATION

	Payload Size	Litres	Outer	Exterior Dimensions WxDxH		Interior Dimensions WxDxH		Tare Weight	
				mm	inches	mm	inches	kg	lbs
Mobile	Series 4 248	2L	White Plastic Box	235 x 210 x 203	9.25 x 8.25 x 8	152 x 127 x 108	6 x 5 x 4.25	3.3	7.2
	Series S OGHC	2L	Standard Bag	235 x 210 x 203	9.25 x 8.25 x 8	152 x 127 x 108	6 x 5 x 4.25	3.5	7.8
	Series 4C 272	2L	White Plastic Box	276 x 251 x 241	10.88 x 9.88 x 9.5	152 x 127 x 108	6 x 5 x 4.25	4.4	9.8
	Series 4 224	2L	Kraft Box	311 x 279 x 286	12.25 x 11 x 11.25	152 x 127 x 108	6 x 5 x 4.25	4.1	9
	Series 4 248	2L	White Plastic Box	276 x 251 x 241	10.88 x 9.88 x 9.5	152 x 127 x 108	6 x 5 x 4.25	4.4	9.6
Medic	Series 4 GHMED	N/A	Standard Bag	356 x 229 x 152	14 x 9 x 6	N/A	N/A	2.7*	6*
One	Series 4 One-V	0.5L	Standard Bag	152 x 127 x 203	6 x 5 x 8	N/A	N/A	0.9	2
	Series 4 One-F	0.5L	Flask	114 x 114 x 247	4.5 x 4.5 x 9.75	N/A	N/A	1	2.2

*Fully loaded container with blood product payload



GOLDEN HOUR™ MOBILE

PCM COOLANT
DURATION 72 HOURS
TEMPERATURE RANGE
2°C TO 8°C
15°C TO 25°C
VOLUME - 2L



GOLDEN HOUR™ MEDIC

PCM COOLANT
DURATION 72 HOURS
TEMPERATURE RANGE
2°C TO 8°C
VOLUME - 2L



GOLDEN HOUR™ ONE - V

PCM COOLANT
DURATION 48 HOURS
TEMPERATURE RANGE
2°C TO 8°C
VOLUME - 0.5L



GOLDEN HOUR™ ONE - F

PCM COOLANT
DURATION 48 HOURS
TEMPERATURE RANGE
2°C TO 8°C
VOLUME - 0.5L

FOR MORE INFORMATION ABOUT PELI BIOTHERMAL PRODUCTS, PLEASE CONTACT:

EMEA OFFICE - INFOEMEA@PELIBIOTHERMAL.COM T: +44 (0)1525 243 770

AMERICAS OFFICE - INFOAMERICAS@PELIBIOTHERMAL.COM T: +1(763) 412-4800 TOLL FREE US: T: +(877) 537-9800

ASIA OFFICE - INFOASIA@PELIBIOTHERMAL.COM T: +65 6681 0095

PELIBIOTHERMAL.COM

ISO17025, ISO9001-2008 Accredited

Golden Hour™, TIC™, Crêdo™, Crêdo Cube™, Crêdo ProEnvision™, Crêdo Xtreme™, Crêdo Duracube™, Crêdo ProMed™, Sherpa Systems™, CoolGuard™ and CoolPall™ are registered trademarks of Pelican Products Inc., its subsidiaries and / or affiliates



©2022 Peli BioThermal Ltd all rights reserved